

Gm130 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q08379

Item Number: CM0002244



Introduction

Rabbit monoclonal antibody targeting human GM130 (GOLGA2, Golgin-95) for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse. Ideal for Golgi apparatus and mitosis research.

[Learn More](#)

Parameter	Value
Product Name	GM130 Rabbit mAb
Remarks/Alias	GM130; DEDHMB
Species	Human
GeneID (Human)	2801
GeneID	2801
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GM130 (Q08379).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q08379
Protein Weight	113kDa
Shipping	Ice bag