

Gna12 Rabbit Polyclonal Antibody (Pab) For Wb, Ihc-P, Elisa - Q03113

Item Number: CM0018075



Introduction

High-quality rabbit polyclonal antibody targeting GNA12 (G protein subunit alpha 12). Suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Recombinant protein immunogen. Ideal for signaling and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	GNA12 Rabbit pAb
Remarks/Alias	RMP; gep; NN3; HG1M1; GNA12
Species	Human
GeneID (Human)	2768
GeneID	2768
Immunogen	Recombinant fusion protein containing amino acids 285-381 of human GNA12 (NP_031379.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03113
Protein Weight	44kDa
Shipping	Ice bag