

Gna14 Rabbit Pab - O95837

Item Number: CM0015851

Introduction

Rabbit polyclonal antibody targeting human GNA14, suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying G protein-coupled signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	GNA14 Rabbit pAb
Remarks/Alias	HG11; GNA14
Species	Human
Gene ID (Human)	9630
Gene ID	9630
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-350 of human GNA14 (NP_004288.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95837
Protein Weight	42kDa
Shipping	Ice bag