


## Gnai2 Rabbit Mab For Wb, Ihc-P, Elisa - P04899

Item Number: CM0002877



### Introduction

Recombinant GNAI2 Rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human and rat GNAI2 (40kDa), a G protein subunit inhibiting adenylate cyclase. Ideal for GPCR signaling and inflammation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | GNAI2 Rabbit mAb                                                                                                      |
| Remarks/Alias    | GIP; HG1C; GNAI2B; H_LUCA15.1; H_LUCA16.1; GNAI2                                                                      |
| Species          | Human                                                                                                                 |
| GeneID (human)   | 2771                                                                                                                  |
| GeneID           | 2771                                                                                                                  |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 250-355 of human GNAI2 (P04899). |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | WB, IHC-P, ELISA                                                                                                      |
| Cross Reactivity | Human, Rat                                                                                                            |
| SWISS            | P04899                                                                                                                |
| Protein Weight   | 40kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |