

Gnas Rabbit Pab For Wb, If/Icc, Ip, Elisa - O95467 / P63092 / P84996 / Q5Jwf2

Item Number: CM0024914

Introduction

Rabbit polyclonal antibody targeting GNAS, suitable for WB, IF/ICC, IP, and ELISA. Reacts with human, mouse, and rat samples. Essential for studies on GNAS-related signaling and neuroendocrine function.

[Learn More](#)

Parameter	Value
Product Name	GNAS Rabbit pAb
Remarks/Alias	AHO; GSA; GSP; POH; GPSA; NESP; SCG6; SgVI; GNAS1; PITA3; C20orf45; GNAS
Species	Human
GenelD (Human)	2778
GenelD	2778
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95467/P63092/P84996/Q5JWF2
Protein Weight	28kDa/44kDa/45kDa/77kDa/109kDa/111kDa
Shipping	Ice bag