


Gne Rabbit Monoclonal Antibody For Wb, Elisa - Q9Y223

Item Number: CM0002923



Introduction

Rabbit monoclonal antibody targeting human GNE (Q9Y223), validated for WB and ELISA. Detects human and rat GNE, a key enzyme in sialic acid biosynthesis, with a molecular weight of 79 kDa.

[Learn More](#)

Parameter	Value
Product Name	GNE Rabbit mAb
Remarks/Alias	NM; DMRV; IBM2; Uae1; GLCNE; GNE
Species	Human
GeneID (human)	10020
GeneID	10020
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 623-722 of human GNE (Q9Y223).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9Y223
Protein Weight	79 kDa
Shipping	Ice bag