

Gnp1 Rabbit Polyclonal Antibody (Oryza Sativa) - Q10Ac7

Item Number: CM0022691

Introduction

GNP1 Rabbit Polyclonal Antibody for Oryza sativa, validated in WB and ELISA. Recombinant immunogen, cross-reacts with Oryza sativa. 42kDa protein. Ideal for plant growth and development research.

[Learn More](#)

Specification	Detail
Product Name	GNP1 Rabbit pAb
Alias	GNP1; OsGA20ox1; qEPD2
Species	Oryza sativa
Gene ID	4334841
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of Oryza sativa GNP1 (Q10AC7)
Source	Rabbit
Category	Polyclonal Antibodies
Applications	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q10AC7
Protein Weight	42kDa
Shipping	Ice bag