


Gnrh1 Rabbit Polyclonal Antibody For Wb And Elisa - P01148

Item Number: CM0026489



Introduction

This rabbit polyclonal antibody against human GnRH1 (also known as LHRH, GRH) is validated for WB and ELISA, cross-reacts with mouse, and is ideal for investigating GnRH-mediated gonadotropin secretion pathways.

[Learn More](#)

Parameter	Value
Product Name	GnRH1 Rabbit pAb
Remarks/Alias	GRH; GNRH; LHRH; LNRH
Species	Human
GeneID (Human)	2796
GeneID	2796
Immunogen	Synthetic peptide corresponding to amino acids 1-92 of human GNRH1 (NP_001076580.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01148
Protein Weight	10kDa
Shipping	Ice bag