

Gns Rabbit Polyclonal Antibody For Wb And Elisa - P15586

Item Number: CM0025365

Introduction

Rabbit polyclonal antibody against human GNS (N-acetylglucosamine-6-sulfatase), validated for WB and ELISA applications. Reacts with human, mouse, and rat. Suitable for glycobiology and lysosomal storage disease research.

[Learn More](#)

Parameter	Value
Product Name	GNS Rabbit pAb
Remarks/Alias	G6S; GNS
Species	Human
Gene ID (Human)	2799
Gene ID	2799
Immunogen	Recombinant protein corresponding to amino acids 293-552 of human GNS (NP_002067.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15586
Protein Weight	62kDa
Shipping	Ice bag