

Gorasp2 Rabbit Pab For Wb, If/Icc, Elisa - Q9H8Y8

Item Number: CM0016478



Introduction

Rabbit polyclonal antibody targeting GORASP2 (GRASP55), a key Golgi reassembly-stacking protein; suitable for WB, IF/ICC, ELISA on human, mouse, rat. Immunogen: recombinant human GORASP2 (aa 353-452).

[Learn More](#)

Parameter	Value
Product Name	GORASP2 Rabbit pAb
Remarks/Alias	p59; GRS2; GOLPH6; GRASP55; GORASP2
Species	Human
GeneID (Human)	26003
GeneID	26003
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 353-452 of human GORASP2 (NP_056345.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8Y8
Protein Weight	47kDa
Shipping	Ice bag