

Gp6 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hcn6

Item Number: CM0023286



Introduction

High-quality GP6 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, rat GP6. Ideal for platelet glycoprotein VI studies in thrombosis and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	GP6 Rabbit pAb
Remarks/Alias	GP1V; GPVI; BDPLT11; GP6
Species	Human
Gene ID (Human)	51206
Gene ID	51206
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 21-267 of human GP6 (NP_057447.5))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HCN6
Protein Weight	37kDa
Shipping	Ice bag