


Gpam Rabbit Polyclonal Antibody For Wb, Elisa - Q9Hcl2

Item Number: CM0026070



Introduction

High-quality rabbit polyclonal antibody against human GPAM (glycerol-3-phosphate acyltransferase 1, mitochondrial). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	GPAM Rabbit pAb
Remarks/Alias	GPAT; GPAT1; GPAM
Species	Human
GeneID (Human)	57678
GeneID	57678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 594-828 of human GPAM (NP_001231878.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9HCL2
Protein Weight	94 kDa
Shipping	Ice bag