

Gprin1 Rabbit Polyclonal Antibody For Wb And Elisa - Q7Z2K8

Item Number: CM0021657

Introduction

High-quality GPRIN1 Rabbit Polyclonal Antibody suitable for western blot and ELISA applications. Reacts with human, mouse, and rat GPRIN1, a protein involved in neurite outgrowth regulation.

[Learn More](#)

Parameter	Value
Product Name	GPRIN1 Rabbit pAb
Remarks/Alias	GRIN1; GPRIN1
Species	Human
GeneID (Human)	114787
GeneID	114787
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human GPRIN1 (NP_443131.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z2K8
Protein Weight	102kDa
Shipping	Ice bag