


Gps1 Rabbit Pab - Q13098

Item Number: CM0025445



Introduction

High-quality rabbit polyclonal antibody against GPS1 (CSN1/COPS1), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for ubiquitin-proteasome pathway and signalosome research. 56kDa.

[Learn More](#)

Parameter	Value
Product Name	GPS1 Rabbit pAb
Remarks/Alias	CSN1; SGN1; COPS1; GPS1
Species	Human
Gene ID (Human)	2873
Gene ID	2873
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 258-527 of human GPS1 (NP_997657.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13098
Protein Weight	56kDa
Shipping	Ice bag