

Granulin Rabbit Polymab® - P28799

Item Number: CM0004096



Introduction

Anti-Granulin Rabbit PolymAb® for WB, IF/ICC, ELISA. Reacts with Human, Mouse. Recombinant monoclonal antibody targeting human Granulin (PGRN, GP88). Ideal for lysosomal function and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Granulin Rabbit PolymAb®
Remarks/Alias	GEP; GP88; PEPI; PGRN; CLN11; PCDGF
Species	Human
GeneID (Human)	2896
GeneID	2896
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 494-593 of human Granulin (NP_002078.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P28799
Protein Weight	64kDa
Shipping	Ice bag