


## Grb10 Rabbit Polyclonal Antibody For Wb, Elisa - Q13322

Item Number: CM0015849



### Introduction

GRB10 rabbit polyclonal antibody validated for WB and ELISA. Detects human GRB10, a key adapter in insulin/IGF1R signaling and Wnt pathways. Ideal for metabolic and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | GRB10 Rabbit pAb                                                                                      |
| Remarks/Alias    | RSS; IRBP; MEG1; GRB-IR; Grb-10; GRB10                                                                |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 2887                                                                                                  |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 1-100 of human GRB10 (NP_005302.3) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, ELISA                                                                                             |
| Cross Reactivity | Human                                                                                                 |
| SWISS            | Q13322                                                                                                |
| Protein Weight   | 67kDa                                                                                                 |
| Shipping         | Ice bag                                                                                               |