

Grb2 Rabbit Monoclonal Antibody For Wb - P62993

Item Number: CM0002347

Introduction

GRB2 Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, rat GRB2. Ideal for studying growth factor receptor signaling, T-cell activation, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	GRB2 Rabbit mAb
Remarks/Alias	ASH; Grb3-3; MST084; NCKAP2; MSTP084; EGFRBP-GRB2; GRB2
Species	Human
GeneID (Human)	2885
GeneID	2885
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 118-217 of human GRB2 (P62993).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62993
Protein Weight	25kDa
Shipping	Ice bag