

Grhl2 Rabbit Polyclonal Antibody For Wb, Elisa - Q6Isb3

Item Number: CM0020177



Introduction

GRHL2 Rabbit Polyclonal Antibody validated for Western blot and ELISA in human and rat. Targets Grainyhead-like protein 2 homolog, a transcription factor crucial for epithelial development and neural tube closure. Ideal for studies in neurulation, epithelial morphogenesis, and keratinocyte biology.

[Learn More](#)

Parameter	Value
Product Name	GRHL2 Rabbit pAb
Remarks/Alias	GRHL2; BOM; DFNA28; ECTDS; TFCP2L3; grainyhead-like protein 2 homolog
Species	Human
GeneID (human)	79977
GeneID	79977
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 110-210 of human GRHL2.
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q6ISB3
Protein Weight	69kDa/71kDa
Shipping	Ice bag