

Grik2 Rabbit Mab For Wb, Ihc-P, Elisa - Q13002

Item Number: CM0012529



Introduction

Recombinant rabbit monoclonal antibody targeting human GRIK2 (GluR6, EAA4). Validated for WB, IHC-P, and ELISA with reactivity across human, mouse, and rat. Ideal for ionotropic glutamate receptor research.

[Learn More](#)

Parameter	Value
Product Name	GRIK2 Rabbit mAb
Remarks/Alias	EAA4; GLR6; MRT6; GLUK6; GLUR6; GluK2; NEDLAS; GRIK2
Species	Human
GeneID (Human)	2898
GeneID	2898
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human GRIK2 (Q13002)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13002
Protein Weight	103kDa
Shipping	Ice bag