


Grp75/Mot/Hspa9 Rabbit Mab For Wb, Elisa - P38646

Item Number: CM0002246



Introduction

High-quality rabbit monoclonal antibody against human Grp75/MOT/HSPA9 (P38646). Validated for WB and ELISA, reacts with human, mouse, rat. Ideal for mitochondrial stress and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	Grp75/MOT/HSPA9 Rabbit mAb
Remarks/Alias	CSA; MOT; MOT2; SAAN; CRP40; EVPLS; GRP75; PBP74; GRP-75; HSPA9B; SIDBA4; MTHSP75; HEL-S-124m; Grp75/MOT/HSPA9
Species	Human
GeneID (Human)	3313
GeneID	3313
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 580-679 of human Grp75/MOT/HSPA9 (P38646)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P38646
Protein Weight	74kDa
Shipping	Ice bag