


Gsk3B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P49841

Item Number: CM0019009



Introduction

GSK3 β rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, rat. Ideal for research on Wnt signaling, glucose metabolism, and Alzheimer's disease.

[Learn More](#)

Parameter	Value
Product Name	GSK3 β Rabbit pAb
Remarks/Alias	GSK3B; gsk-3 β ; GSK3 β
Species	Human
Gene ID (Human)	2932
Gene ID	2932
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49841
Protein Weight	47kDa
Shipping	Ice bag