

Gsr Rabbit Polyclonal Antibody For Wb And Elisa - P00390

Item Number: CM0012886



Introduction

Rabbit polyclonal antibody against human glutathione reductase (GSR) validated for WB and ELISA, cross-reactive with mouse. Suitable for studying redox biology, oxidative stress, and hemolytic anemia pathways.

[Learn More](#)

Parameter	Value
Product Name	GSR Rabbit pAb
Remarks/Alias	GR; GSRD; HEL-75; HEL-S-122m; GSR
Species	Human
GeneID (Human)	2936
GeneID	2936
Immunogen	Recombinant protein corresponding to amino acids 273-522 of human GSR (NP_000628.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00390
Protein Weight	56kDa
Shipping	Ice bag