


Gsta4 Rabbit Mab For Wb, Elisa - O15217

Item Number: CM0005980



Introduction

Recombinant rabbit monoclonal antibody targeting human GSTA4 (GTA4, GSTA4-4), validated for Western blot (WB) and ELISA. Detects a 26kDa protein; cross-reacts with human samples. Ideal for studying glutathione S-transferase A4 function and 4-HNE detoxification.

[Learn More](#)

Parameter	Value
Product Name	GSTA4 Rabbit mAb
Remarks/Alias	GTA4; GSTA4-4; GSTA4
Species	Human
GeneID (human)	2941
GeneID	2941
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15217
Protein Weight	26kDa
Shipping	Ice bag