


Gstt1 Rabbit Monoclonal Antibody For Wb, Elisa - P30711

Item Number: CM0004219



Introduction

GSTT1 rabbit monoclonal antibody validated for Western blot (WB) and ELISA. Reacts with human, mouse, rat. Immunogen: recombinant GSTT1 (aa 1-240). For glutathione S-transferase studies.

[Learn More](#)

Parameter	Value
Product Name	GSTT1 Rabbit mAb
Remarks/Alias	GSTT1
Species	Human
GeneID	2952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human GSTT1. (NP_000844.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30711
Protein Weight	27kDa
Shipping	Ice bag