


Gtf2F1 Rabbit Polyclonal Antibody For Wb, Elisa - P35269

Item Number: CM0020789



Introduction

GTF2F1 Rabbit Polyclonal Antibody (P35269) validated for WB and ELISA. Recognizes human, mouse, rat. General transcription factor IIF subunit 1 (TFIIF-alpha), 58kDa, rabbit source.

[Learn More](#)

Parameter	Value
Product Name	GTF2F1 Rabbit pAb
Remarks/Alias	BTf4; RAP74; TF2F1; TFIIF; GTF2F1
Species	Human
Gene ID (Human)	2962
Gene ID	2962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 390-517 of GTF2F1 (NP_002087.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35269
Protein Weight	58kDa
Shipping	Ice bag