


GucY2C Rabbit Polyclonal Antibody For Wb And Elisa - P25092

Item Number: CM0017968



Introduction

Rabbit polyclonal antibody targeting human GUCY2C (guanylyl cyclase C), validated for Western blot and ELISA applications. Cross-reacts with mouse. Ideal for enterotoxin signaling and gastrointestinal disease research.

[Learn More](#)

Parameter	Value
Product Name	GUCY2C Rabbit pAb
Remarks/Alias	GC-C; STAR; DIAR6; GUC2C; MECIL; MUCIL; GUCY2C
Species	Human
Gene ID (Human)	2984
Gene ID	2984
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-330 of human GUCY2C (NP_004954.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25092
Protein Weight	123kDa
Shipping	Ice bag