

Habp2 Rabbit Pab For Wb, If/Icc, Elisa - Q14520

Item Number: CM0030678

Introduction

Rabbit polyclonal antibody targeting HABP2 (Factor VII-activating protease) for applications including WB, IF/ICC, and ELISA. Detects human and mouse HABP2. Suitable for coagulation and tumor suppression studies.

[Learn More](#)

Parameter	Value
Product Name	HABP2 Rabbit pAb
Remarks/Alias	FSAP; HABP; PHBP; HGFAL; NMTC5; HABP2
Species	Human
GeneID (human)	3026
GeneID	3026
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-313 of human HABP2 (NP_004123.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14520
Protein Weight	63kDa
Shipping	Ice bag