

## Haus8 Rabbit Pab For Wb - Q9Bt25

Item Number: CM0027364



### Introduction

Validated for WB, IF/ICC, and ELISA, this HAUS8 rabbit polyclonal antibody recognizes human and mouse HAUS8 (Q9BT25), supporting research on mitotic spindle assembly, centrosome integrity, and cytokinesis.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | HAUS8 Rabbit pAb                                                                                                    |
| Remarks/Alias    | DGT4; HICE1; NY-SAR-48; HAUS8                                                                                       |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 93323                                                                                                               |
| Gene ID          | 93323                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 166-410 of human HAUS8 (NP_219485.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse                                                                                                        |
| SWISS            | Q9BT25                                                                                                              |
| Protein Weight   | 45kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |