

Havcr1 Rabbit Pab - Q96D42

Item Number: CM0030266



Introduction

Rabbit polyclonal antibody against human HAVCR1 (TIM-1/KIM-1) for WB and ELISA, reacting with human, mouse, and rat. Useful for studies on kidney injury, viral entry, and T-cell immunity.

[Learn More](#)

Parameter	Value
Product Name	HAVCR1 Rabbit pAb
Remarks/Alias	TIM; KIM1; TIM1; CD365; HAVCR; KIM-1; TIM-1; TIMD1; TIMD-1; HAVCR-1
Species	Human
GeneID (Human)	26762
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-125 of human HAVCR1 (NP_036338.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96D42
Protein Weight	39kDa
Shipping	Ice bag