


Hbg1/2 Rabbit Mab (Wb, Ihc, If, Elisa) - P69891/P69892

Item Number: CM0003094



Introduction

High-specificity rabbit monoclonal antibody against human HBG1/2. Validated for WB, IHC, IF-P, ELISA, detects ~12 kDa fetal hemoglobin gamma chain. Ideal for developmental hematology and HbF research.

[Learn More](#)

Parameter	Value
Product Name	HBG1/2 Rabbit mAb
Remarks/Alias	HBG-T2; HBGA; HBGR; HSGGL1; PRO2979; HBG1/2
Species	Human
GeneID	3047/3048
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HBG1 (P69891).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P69891/P69892
Protein Weight	12 kDa
Shipping	Ice bag