


Hbg1/Hbg2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P69892

Item Number: CM0019407



Introduction

HBG1/HBG2 Rabbit Polyclonal Antibody detects fetal hemoglobin gamma chains in WB, IF/ICC, ELISA. Ideal for hemoglobinopathy research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	HBG1/HBG2 Rabbit pAb
Remarks/Alias	TNCY; HBG-T1; HBG1/HBG2
Species	Human
GeneID (Human)	3048
GeneID	3048
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HBG1/HBG2 (NP_000175.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P69892
Protein Weight	16kDa
Shipping	Ice bag