

Hcar3 Rabbit Polyclonal Antibody For Wb - P49019

Item Number: CM0023550



Introduction

High-quality HCAR3 rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying hydroxycarboxylic acid receptor 3 function in adipocyte biology and metabolism.

[Learn More](#)

Parameter	Value
Product Name	HCAR3 Rabbit pAb
Remarks/Alias	HCA3; HM74; PUMAG; Puma-g; GPR109B; HCAR3
Species	Human
GeneID (Human)	8843
GeneID	8843
Immunogen	Synthetic peptide corresponding to amino acids 300-392 of human HCAR3 (NP_006009.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49019
Protein Weight	44kDa
Shipping	Ice bag