

## Hccs Rabbit Pab - P53701

Item Number: CM0029959

### Introduction

High-quality HCCS Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Targets holocytochrome c-type synthase, essential for cytochrome c maturation. Ideal for mitochondrial disorder research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | HCCS Rabbit pAb                                                                                                 |
| Remarks/Alias    | MLS, CCHL, MCOPS7, LSDMCA1, HCCS                                                                                |
| Species          | Human                                                                                                           |
| GeneID (Human)   | 3052                                                                                                            |
| GeneID           | 3052                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human HCCS (NP_005324.3) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IF/ICC, ELISA                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P53701                                                                                                          |
| Protein Weight   | 31kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |