

Hcg Beta Rabbit Monoclonal Antibody For Elisa, Wb, Ihc-P - P0Dn86

Item Number: CM0006686



Introduction

Rabbit monoclonal hCG Beta antibody for ELISA, WB, IHC-P. Specific to human chorionic gonadotropin beta subunit. Recombinant immunogen; 18kDa. Essential for pregnancy research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	hCG Beta Rabbit mAb
Remarks/Alias	CGB; LHB; CGB5; CGB7; CGB8; hCGB; hCG Beta
Species	Human
GeneID (Human)	1082
GeneID	1082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-165 of human hCG Beta1 (NP_000728.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	ELISA, WB, IHC-P
Cross Reactivity	Human
SWISS	P0DN86
Protein Weight	18kDa
Shipping	Ice bag