


Hcn1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O60741

Item Number: CM0021537



Introduction

High-quality HCN1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Targets human and mouse potassium/sodium hyperpolarization-activated cyclic nucleotide-gated channel 1, involved in pacemaker currents and neuronal excitability.

[Learn More](#)

Parameter	Value
Product Name	HCN1 Rabbit pAb
Remarks/Alias	BCNG1; DEE24; HAC-2; BCNG-1; EIEE24; GEFSP10; HCN1
Species	Human
Gene ID	348980
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 791-890 of human HCN1 (NP_066550.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O60741
Protein Weight	99kDa
Shipping	Ice bag