

Hdac1 Mouse Mab - Q13547

Item Number: CM0003280

Introduction

Monoclonal mouse antibody targeting human HDAC1 (UniProt Q13547). Validated for Western blot and ELISA applications. Suitable for research, diagnostic assay development, and IVD raw material procurement.

[Learn More](#)

Parameter	Value
Product Name	HDAC1 Mouse mAb
Gene ID (Human)	3065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 393-482 of human HDAC1 (Q13547)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS Accession	Q13547
Protein Weight	62kDa
Shipping	Ice bag