


## Hdac2 Rabbit Pab - Q92769

Item Number: CM0019668



### Introduction

Rabbit polyclonal antibody against HDAC2 (histone deacetylase 2) for WB, IHC-P, IP, ELISA. Detects human, mouse, rat HDAC2 (~55kDa). Suitable for epigenetic and chromatin research.

[Learn More](#)

| Parameter        | Value                                                                                                 |
|------------------|-------------------------------------------------------------------------------------------------------|
| Product Name     | HDAC2 Rabbit pAb                                                                                      |
| Remarks/Alias    | HD2; RPD3; YAF1; KDAC2; HDAC2                                                                         |
| Species          | Human                                                                                                 |
| GeneID (Human)   | 3066                                                                                                  |
| GeneID           | 3066                                                                                                  |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 390-488 of human HDAC2 (NP_001518.3) |
| Source           | Rabbit                                                                                                |
| Category         | Polyclonal Antibodies                                                                                 |
| Application      | WB, IHC-P, IP, ELISA                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                     |
| SWISS            | Q92769                                                                                                |
| Protein Weight   | 55kDa                                                                                                 |
| Shipping         | Ice bag                                                                                               |