

Hdac4 Rabbit Mab - P56524

Item Number: CM0002195



Introduction

High-quality HDAC4 rabbit monoclonal antibody validated for WB and ELISA. Specific for human HDAC4, a 119kDa class II histone deacetylase essential in epigenetic regulation, muscle maturation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HDAC4 Rabbit mAb
Remarks/Alias	HD4; AHO3; BDMR; HDACA; HA6116; HDAC-4; HDAC-A; NEDCHF; NEDCHID; HDAC4
Species	Human
GeneID (Human)	9759
GeneID	9759
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HDAC4 (P56524).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P56524
Protein Weight	119kDa
Shipping	Ice bag