


Hdac5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uql6

Item Number: CM0027130



Introduction

Rabbit polyclonal antibody targeting human HDAC5, a histone deacetylase involved in epigenetic regulation and muscle development. Validated for WB, IF/ICC, and ELISA. Ideal for transcriptional repression studies.

[Learn More](#)

Parameter	Value
Product Name	HDAC5 Rabbit pAb
Remarks/Alias	HD5; NY-CO-9; HDAC5
Species	Human
Gene ID (Human)	10014
Gene ID	10014
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-450 of human HDAC5 (NP_005465.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9UQL6
Protein Weight	122kDa
Shipping	Ice bag