

Hdac6 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q9Z2V5

Item Number: CM0027894

Introduction

Rabbit polyclonal anti-HDAC6 antibody for WB, IF/ICC, IHC-P, and ELISA. Detects human, mouse, and rat HDAC6. Ideal for research on tubulin deacetylation, autophagy, and aggresome pathways.

[Learn More](#)

Parameter	Value
Product Name	HDAC6 Rabbit pAb
Remarks/Alias	Hd6; Sfc6; Hdac5; mHDA2
Species	Mouse
GeneID (Human)	10013
GeneID	15185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 836-1037 of mouse HDAC6 (NP_001123888.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Z2V5
Protein Weight	126kDa
Shipping	Ice bag