


Hdac6 Rabbit Polyclonal Antibody (Wb, Ihc-P, Elisa) - Q9Ubn7

Item Number: CM0021562



Introduction

HDAC6 Rabbit Polyclonal Antibody validated in WB, IHC-P, ELISA. Detects human, mouse, rat HDAC6 (131 kDa). Suitable for studies on tubulin deacetylation, autophagy, and cancer.

[Learn More](#)

Parameter	Value
Product Name	HDAC6 Rabbit pAb
Remarks/Alias	HD6; JM21; CPBHM; PPP1R90; HDAC6
Species	Human
GeneID	10013
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HDAC6 (NP_006035.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBN7
Protein Weight	131kDa
Shipping	Ice bag