


Hdac8 Rabbit Mab For Ihc-P, If/Icc, Ip, Elisa - Q9By41

Item Number: CM0008366



Introduction

HDAC8 Rabbit mAb validated for IHC-P, IF/ICC, IP, ELISA with human, mouse, rat reactivity. Targets histone deacetylase 8, a key chromatin regulator involved in transcription and cell cycle. Suitable for epigenetics and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HDAC8 Rabbit mAb
Remarks/Alias	HD8; WTS; RPD3; CDA07; CDLS5; KDAC8; MRXS6; HDACL1; C8
Species	Human
GeneID (Human)	55869
GeneID	55869
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HDAC8 (Q9BY41).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY41
Protein Weight	42kDa
Shipping	Ice bag