


Hdac9 Rabbit Mab - Q9Ukv0

Item Number: CM0006510



Introduction

High-quality HDAC9 Rabbit monoclonal antibody for Western blot, IHC-P, IF/ICC, ELISA. Detects human and rat HDAC9. Ideal for epigenetic and transcriptional regulation research - Q9UKV0.

[Learn More](#)

Parameter	Value
Product Name	HDAC9 Rabbit mAb
Remarks/Alias	HD7; HD9; HD7b; HDAC; HDRP; MITR; HDAC7; HDAC7B; HDAC9B; HDAC9FL; HDAC9
Species	Human
GeneID (Human)	9734
GeneID	9734
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UKV0
Protein Weight	111kDa
Shipping	Ice bag