

Hdgfrp2 Rabbit Polyclonal Antibody - Hdgfrp2

Item Number: CM0025514

Introduction

Rabbit polyclonal antibody targeting human HDGFRP2, validated in WB, IF/ICC, and ELISA. Cross-reacts with mouse. Immunogen corresponds to C-terminal region (aa 501-670). Ideal for research on hepatoma-derived growth factor-related protein 2.

[Learn More](#)

Parameter	Value
Product Name	HDGFRP2 Rabbit pAb
Remarks/Alias	HRP2; HDGF2; HRP-2; HDGF-2; HDGFRP2
Species	Human
Gene ID	84717
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 501-670 of human HDGFRP2 (NP_116020.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
Protein Weight	74kDa
Shipping	Ice bag