

Hdgfrp3 Rabbit Polyclonal Antibody - Hdgfrp3

Item Number: CM0026730



Introduction

Rabbit polyclonal antibody targeting hepatoma-derived growth factor-related protein 3 (HDGFRP3/HDGF2/HRP-3). Validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for cancer and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	HDGFRP3 Rabbit pAb
Remarks/Alias	HDGF2; HRP-3; HDGF-2; CGI-142; HDGFRP3
Species	Human
Gene ID	50810
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-203 of human HDGFRP3 (NP_057157.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	22kDa
Shipping	Ice bag