

Heme Oxygenase 1 (Ho-1/Hmox1) Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P09601

Item Number: CM0013872



Introduction

Rabbit polyclonal antibody targeting Heme Oxygenase 1 (HO-1/HMOX1), validated for WB, IF-P, IHC-P, and ELISA in human, mouse, and rat samples. Suitable for heme metabolism, cancer, and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	Heme Oxygenase 1 (HO-1/HMOX1) Rabbit pAb
Remarks/Alias	HO-1; HSP32; HMOX1D; bK286B10; Heme Oxygenase 1 (HO-1/HMOX1)
Species	Human
GeneID (Human)	3162
GeneID	3162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human Heme Oxygenase 1 (HO-1/HMOX1) (NP_002124.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09601
Protein Weight	33kDa
Shipping	Ice bag