

Hes5 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q5Ta89

Item Number: CM0011972

Introduction

HES5 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat HES5. Ideal for neurogenesis and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	HES5 Rabbit mAb
Remarks/Alias	bHLHb38; HES5
Species	Human
GeneID (Human)	388585
GeneID	388585
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 12-142 of human HES5 (Q5TA89)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5TA89
Protein Weight	18kDa
Shipping	Ice bag