


Hexa Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P06865

Item Number: CM0030222



Introduction

HEXA rabbit polyclonal antibody targeting beta-hexosaminidase subunit alpha. Suitable for WB, IF-P, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for lysosomal storage disease and lipid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	HEXA Rabbit pAb
Remarks/Alias	TSD; HEXA
Species	Human
Gene ID (Human)	3073
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 270-529 of human HEXA (NP_000511.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06865
Protein Weight	61kDa
Shipping	Ice bag