


Hfe2 Rabbit Pab - Q6Zvn8

Item Number: CM0029504



Introduction

HFE2 Rabbit pAb is a polyclonal antibody validated for Western blot and ELISA, targeting human hemojuvelin with cross-reactivity to mouse and rat, for iron homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	HFE2 Rabbit pAb
Remarks/Alias	JH; HFE2; RGMG; HFE2A
Species	Human
GeneID	148738
Immunogen	Recombinant protein corresponding to amino acids 1-174 of human HFE2 (NP_998817.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6ZVN8
Protein Weight	45kDa
Shipping	Ice bag