


Hibch Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q6Nvy1

Item Number: CM0012035



Introduction

HIBCH Rabbit Monoclonal Antibody validated for WB, IHC-P, ELISA in Human, Mouse, Rat. Detects 43kDa mitochondrial 3-hydroxyisobutyryl-CoA hydrolase; ideal for valine catabolism and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	HIBCH Rabbit mAb
Remarks/Alias	HIBYLCOAH
Species	Human
Gene ID (Human)	26275
Gene ID	26275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-330 of human HIBCH (NP_055177.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6NVY1
Protein Weight	43kDa
Shipping	Ice bag